



The internal diameter needs to match the retention zone (red line)--or area on the male post the O-Ring rests--to assure proper fit.

[OSO O-Ring Instructions](#)

[OSO O-Ring Components](#)

OSO Replacement Rings Size Chart

OSO Ring	Sizing	Internal Diameter (ID)		Outside Diameter (OD)		Price
		inch	mm	inch	mm	
White #2 Pkg 12	Mini	3/64"	1.20 mm	9/64"	3.57 mm	\$21.00
White #3 Pkg 12	Regular	4/64"	1.60 mm	12/64"	4.78mm	\$21.00
Black #3 Pkg 6	Regular	4/64"	1.60 mm	12/64"	4.78mm	\$24.00
Black #4 Pkg 6	Large	5/64"	1.98 mm	13/64"	5.16 mm	\$27.00

Black #5 Pkg 6	Extra Large	6/64"	2.39 mm	14/64"	5.55 mm	\$27.00
			Metal Retainer Rings #2 & #3 (ea)			\$21.00
			Metal Retainer Rings #4 (ea)			\$24.00
			Metal Retainer Rings #5 (ea)			\$25.00